PORTABLE REAL-TIME PCR FOR RAPID SCREENING

SHINEWAY TECHNOLOGY CORPORATE



The Company

Shineway Technology Corporate was founded by a team of researchers at Hong Kong University of Science and Technology in September 2016.

The company is a national high-tech enterprise specializing in the development of **rapid nucleic acid molecular in vitro diagnostic technology**, which is committed to research and development, production and sales of world-class on-site rapid molecular in vitro diagnostic medical device products.



COVID-19 CORONAVIRUS

Coronavirus cases soar to > 70,000.

China reported that 1,716 medical workers were infected. (1,502 of those cases in Hubei Province), and six of these workers have died, the Times reports.

Hubei Province, where the outbreak started, has changed how it diagnoses the virus, causing the number of cases to skyrocket in the area.

About 74,456 confirmed coronavirus cases (primarily in mainland China), according to the Johns Hopkins virus dashboard.

2,800 deaths due to the infection. Number of deaths worldwide is leading over SARS.

The coronavirus cruise ship quarantine is a scary public health experiment

The second-largest outbreak outside China is floating at sea — and health experts are increasingly worried.

By Julia Belluz | @juliaoftoronto | julia.belluz@voxmedia.com | Updated Feb 12, 2020, 8:33am EST





Their concern is driven by the details emerging from the ship, which paint an alarming portrait. According to a **must-read New York Times** report, as of February 10, Japanese officials had only tested 439 of the 3,700 passengers onboard for the coronavirus, citing shortages of testing supplies. This means that at least a fraction of the remaining thousands onboard may have the virus and not know it.

There has also been criticism of the slow pace of testing, although officials have said they hope to increase the number of daily tests from 300 to 1,000.

"On the shin infections are getting very dense ... it now provides a

CASE STUDY



NEWS:

About 3,700 people are being quarantined on the Diamond Princess cruise ship off Japan due to COVID-19. 213 were tested positive (as of 16th of Feb 2020), including a onboard Japanese health official.



Problems Statement:

Epidemic outbreak Samples-in queue for testing Lacking of portable laboratory Complicated sample handling Long testing time Long quarantine time




VISION

Remote Laboratory




To provide testing solution without boundaries



Point of Care Testing (POCT) Solution





DEFINITION:

POCT is defined as an investigation taken at the time of the consultation with instant availability of results to make immediate and informed decisions about patient care..."


Kamlesh Khunti, a Professor of Primary Care Diabetes and Vascular Medicine at the University of Leicester

POCT Type Difference



Molecular Test

- Test based on DNA or RNA. High technological barriers
- High specificity and sensitivity
- ✓ Short window period




Immunoassays

- Test based on protein, metabolites, carbohydrate, lipids, Inorganic ion, etc.
- Mature and prevalent
 Low sensitivity and limited range
 Long window period

Viruses can be tested at the primary stage by Molecular Test, which can effectively block the spread of the virus.

Challenges in Diagnosis Industry







Key Product

Microfluidic PCR Reagents

-Nucleic Acid amplification Reaction Components -Specification: 24T/kit

Dimension (mm)	330 x 250 x 180
Weight (kg)	2.5


Portable Real-time PCR

Microfluidic PCR Chip

-PCR reaction chamber -Chamber Volume 2-18microlitre

PCR Nucleic Acid Analyzer


-consist of temperature control systems, electro-optical system, microfluidic chip, data processing and overall control chip, and etc.



Supporting jig and fixture



Rapid Detection Process




Test Samples:

Throat Swabs, Nasal Swabs, Deep cough sputum, blood specimens, fecal specimens, conjunctival swab specimens Scenario Application: E.g.: Test Kit for COVID-19: Rapid Screening, total estimated time: 50 mins Pool test or accuracy check with 3rd parties' solution, total estimated time: 1.5 hours


* Note: Total estimated times varies for different test kits with different protocol



Different Protocol for COVID-19



* Note: Total estimated times varies according to the third parties suggested procedure





Portable Real-time PCR

epidemic

Affordable, Convenient and Precise.

Affordable Half the price of the conventional analyzer.

1

2



On-site Rapid PCR Diagnostic System

Move the central laboratory to the front line of the epidemic

- User friendly interface
- Portable for outdoor testing, zero delivery time
- High reaction efficiency, rapid screening
- Less reagent consumption, lower testing cost
- High detection sensitivity
- New light path and image segmentation algorithm design with 99% accuracy
- Support wireless network for easy data storage and transmission

On-site Rapid PCR Diagnostic System

Broad Application Range

- Pathogen in Epidemic Area (COVID-19)
- Flu Virus
- Enterovirus
- **Respiratory Virus**
- Blood Infectious disease (Hep-B, HIV, HPV, etc)

As reference by **Nature Biotechnology** (February 19, 2020) to be the reliable diagnostic

https://www.nature.com/articles/d41587-020-00002-2

Hong Kong University of Science and Technology

On-site rapid molecular diagnostic system based on Shenzhen Shineway Technology

Integrated microfluidic PCR test employing silicon-based micro-heater module for rapid heating and processing of test samples; delivers test results in 40 minutes

Platform has European Union CE approval

Application Area





Non-Clinical -CDC -Custom -Entry-Exit Inspection and Quarantine Bureaus -Inspection and Quarantine -Animal Clinic -Public Safety (Airport, Train Station)



Clinical -Emergency Unit -Hospital -Blood Centre -Clinic




Domestic

-Self Health Check -Enterprise -School -Construction Units -Industrial Park



Product Performance





E. coli test result by Shineway SWM-01


E. coli test result by LightCycler[®] 480

99% Accuracy

Product Performance-COVID-19



99% Accuracy



Fluorescence vs Cycle Graph

COVID-19 test result by Shineway SWM-01

Sample Amount in the chip well of SWM-01: 20microL

Sample Amount in the tube-based by LightCycler[®] 480: 200microL

Certification and Qualification





Our Team Leader





http://www.phys.ust.hk/phwen/

PhD in Applied Physics Chinese Academy of Sciences, 1995.

> Professor Department of Physics, HKUST

Director of Shenzhen Shineway Technology Corporation

DR. WEIJIA WEN

Main Research Interest:


Soft condensed matter physics, electrorheological (ER) and magnetorheological (MR) fluids, field-induced pattern and structure transitions, micro- and nano-fluidic controlling, microsphere and nanoparticle fabrications, thin film physics, band gap materials, metamaterials and nonlinear optical materials.

- ✔ Published more than 300 SCI papers, others cited more than 10,000
- ✓ Obtained 18 US patents, more than 40 Chinese patents
- Hosted more than 40 national, provincial and ministerial levels, innovation funds and other research projects, with a total research funding of more than 80 million yuan


2000 National Outstanding Youth Fund (Class B)

- 2013 Chief scientist of scientific research in the non-profit industry (agriculture)
- 2014 National Natural Science Second Prize
- 2016 Shenzhen Peacock Venture Team Leader

Clients and Collaborations



Clients and Collaboration




Product Comparison with Competitors

	Accuracy	Portability	Fastest Diagnostic Time	Product on Market	Technology Used and Future Potential
Shineway Hong Kong	>99%	 2kg, on-site embedded system no PC required 	40 minutes	 IVD CE was obtained in 2019 Product currently on sale 	 Microfluidic chip PCR technology Thin film micro heater Great future potential
Novacyt UK	99%	4kg, on-sitePC required	<2 hour	 CE was obtained in 2020 For research use only Product launch estimated last quarter of 2020 	 Tube based PCR, Traditional electronic heater Saturated technology
Biomeme USA	70%	 1kg, hand-held Mobile phone/PC required 	1 hour	 Currently prototype build, no further progress 	 Traditional PCR technology Traditional electronic heater More on software development

Roadmap





NEWS UPDATE

China Daily - HKUST's invention to speed up diagnoses of novel coronavirus http://www.chinadailyhk.com/article/120331#HKUST%E2%80%99s-invention-to-speed-up-diagnoses-of-novelcoronavirus

Shenzhen Daily - Coronavirus fast-detection device developed http://www.szdaily.com/content/2020-02/07/content_22836688.htm

Dimsum Daily - Research team from HKUST launches world's fastest portable detector for novel coronavirus https://www.dimsumdaily.hk/research-team-from-hkust-launches-worlds-fastest-portable-detector-for-novelcoronavirus/

Mirage News - HKUST Research Team Invents World's Fastest Coronavirus Detection Device Offering Diagnostic Results https://www.miragenews.com/hkust-research-team-invents-world-s-fastest-coronavirus-detection-device-offeringdiagnostic-results/

Unwire.HK:科大全球最快40分鐘肺炎病毒檢測機 溫維佳教授:售價平一半+政府正商討引入 https://unwire.hk/2020/02/07/cornoavirusdetection/lifetech/health/?fbclid=IwAR3cutVuTmQJuCBVSoile41bHLZebtm PilOpdSXJJ dOguNfXTTNoAoSAQ

文匯報-抗擊新型肺炎:科大手提箱測冠毒快絕全球 http://paper.wenweipo.com/2020/02/07/HS2002070035.htm



THANK YOU

